

PORTABLE X-RAY

eagle. X-RAY

DABI ATLANTE

YOUR EFFICIENCY WITH THIS MOBILITY

ULTRA LOW DOSE TECHNOLOGY



MOBILITY


A product developed to bring efficiency and sophistication to everywhere.

LIGHT AND COMPACT

Its lightness makes it easy to handle, bringing practicality that increases daily productivity.

PERFECT PERFORMANCE

Ready for the use of intraoral sensor, phosphorus plate and radiographic film.



- The Portable X-Ray ensures safety to patient and operator, with mechanisms that reduce the amount of patient exposure.
- The Portable X-Ray minimizes the leakage radiation through a radiopaque material collimator, avoiding unnecessary user exposure and directing the radiation to the exam target.
- Nominal voltage 60kVp, tube current 2.5 mA, and focal spot 0.4x0.4mm.

Illustrative Photos

Illustrative Photos. Dabi Atlante reserves the right to make changes to the design and technical specifications of the products upon authorization by ANVISA to continue improving them. Portable Dental X-Rays: Anvisa: 10101130092. 51000000180: Version - Out/2025.



dabiatlante



dabiatlante_oficial



dabiatlante

DABI ATLANTE