

PORTABLE X-RAY

# eagle. X-RAY

DABI ATLANTE

EFFICIENCY WITH MOBILITY

ULTRA LOW DOSE TECHNOLOGY



## MOBILITY

A product developed to bring efficiency and sophistication - anywhere.

## LIGHT AND COMPACT

Its lightness makes it easy to handle, bringing practicality that increases daily productivity.

## PERFECT PERFORMANCE

Ready for use with an intraoral sensor, phosphorus plate, or radiographic film.



PATIENT PROFILE

ADJUSTMENTS KEY

ON / OFF

X-RAY MODE

Digital Intraoral Sensor  
Phosphorus Plate  
Radiographic Film

TOOTH TYPE



- The Portable X-Ray ensures safety for the patient and operator, with mechanisms that reduce patient exposure.
- The Portable X-Ray minimizes radiation leakage through a radiopaque material collimator, avoiding unnecessary user exposure and directing the radiation to the exam target.
- Nominal voltage: 60kVp  
Tube current: 2.5 mA  
Focal spot: 0.7x0.7mm

Illustrative Photos

Illustrative Photos. Dabi Atlante reserves the right to make changes to the design and technical specifications of the products, upon authorization by ANVISA, to continue improving them. Portable Dental X-Rays: Anvisa: 10101130092. 51000000180: Version - May2023



dabiatlante



dabiatlante\_oficial



dabiatlante

DABI ATLANTE

Dabi Atlante is a member of the Alliage group. | Contact: sales@alliage-international.com **alliage**