


DENTAL UNIT

Personal

DABI ATLANTE TECHNOLOGY **air**

- BI-ARTICULATED HEADREST
- FOLDING ARMREST
- MULTIFUNCTION JOYSTICK PEDAL
- DELIVERY UNIT ARM WITH PNEUMATIC LOCK
- WATER UNIT WITH PROXIMITY SENSOR AND CUP HOLDER
- 5 LED DENTAL LIGHT WITH PROXIMITY SENSOR (OPTIONAL)
- RETRACTABLE RODS WITH LOCKING SYSTEM (HASTEFLX VERSION)



WATER UNIT WITH PROXIMITY SENSOR AND CUP HOLDER

DENTAL UNIT

Personal

hasteflex



DENTAL CHAIR

- 90° Foldable patient support arm
- Anatomic and bi-articulated headrest
- Joystick pedal with return to zero position, 3 working positions, spit position, dental light activation, and movement block.
- Ambidextrous

DELIVERY UNIT

- Arm with pneumatic lock
- Up to 5 terminals
- Retractable rods with locking system (Hasteflex Version)
- Removable stainless steel top (Air Version)
- Composition:
 - 01 triple syringe
 - 01 terminal with spray for low rotation
 - 01 terminal for turbine

WATER UNIT

- Debris filter located at the base of the chair
- 90° rotatable water unit
- Electric water activation
- Venturi suction terminal
- Proximity sensor
- Cup holder

5 LED DENTAL LIGHT WITH PROXIMITY SENSOR (OPTIONAL)

- Lower power consumption and longer durability
- Proximity sensor
- Adjustable light intensity: up to 35.000 LUX
- Handles: Removable and autoclavable
- Frontal Protection: Removable
- Standard Light: 3 LED, with proximity sensor

DABI COLOURS

UPHOLSTERY THAT INSPIRES YOUR WORK ROUTINE

PVC OR VISCOELASTIC*



LEATHER*



*Optional.

OPTIONS

Personal

YOUR DENTAL OFFICE
DESERVES EVOLUTION

N: NOT APPLICABLE | O: OPTIONAL | S: STANDARD

DENTAL CHAIR	PERSONAL AIR	PERSONAL HASTEFLEX
Column Lift System	N	N
Trendelenburg Position	N	N
Massager	N	N
Standard PVC	S	S
Viscoelastic PVC	O	O
Leather	O	O
Bi-articulated Headrest	S	S
Pneumatic Multi-articulated Headrest	N	N
1x Fixed Armrest + 1x 90° Folding Armrest	S	S
1x Fixed Armrest + 1x Pneumatic retractable Armrest	N	N
1x Fixed Armrest + 1x Retractable Armrest (V Version)	N	N
Foldable Leg Support (at knee height)	N	N
Feet Protector	O	O
Foot Control (3 keys)	N	N
Foot Control (5 keys)	S	S
Foot Control (7 keys)	N	N
Foot Control Pedal - with Chip Blower function	N	N
Soft Comfort	O	O
Single Connection Box	O	O
Ambidextrous Dental Chair	O	O

DELIVERY UNIT	PERSONAL AIR	PERSONAL HASTEFLEX
Arm with Air Lock	N	S
Stainless Steel Top	S	S
Upper Tray (Acrylic)	N	N
Lower Tray (Injected)	N	N
Control Pad Kit	O	O
X-ray Viewer	O	O
Flush Kit	O	O
Manometer	N	N
High Speed Terminal	S	S
Low Rotation Terminal	S	S
Low Rotation Spray Terminal	N	N
Triple Syringe	S	S
Triple Injected Syringe	N	N
Syringe Heater Kit	O	O
Number of Terminals in the Equipment	UP TO 5 (3S & 2O)	UP TO 5 (3S & 2O)
Electric Micromotor Kit	O	O
Borden Terminal Kit - with Spray	O	O
Fibre Optic Terminal Kit	O	O
Curing Light Kit	O	O
Ultrasonic Scaler Kit	O	O
Bicarbonate Jet Kit	O	O

WATER UNIT	PERSONAL AIR	PERSONAL HASTEFLEX
Proximity Sensor Kit	O	O
Cup Holder Kit	O	O
Electrical activation of water	N	N
Assistant Arm - without control panel	O	O
Assistant Arm - with control panel	O	O
Venturi Suction Terminal	S	S
Vac Plus Suction Terminal	O	O
Vacuum Pump Suction Terminal	O	O
Triple Injected Syringe	O	O
Triple Syringe Heating	N	N
Water Unit Syringe Heater Kit	O	O

DENTAL LIGHT	PERSONAL AIR	PERSONAL HASTEFLEX
Persus LED	O	O
Sirius Sensor 3 LEDs	S	S
Sirius Sensor 5 LEDs	O	O
Sirius Sensor 5 (3x2) LEDs	O	O

OTHERS	PERSONAL AIR	PERSONAL HASTEFLEX
Multimedia Kit Support	O	O

Illustrative Photos. Colours may have altered in the graphic reproduction. Dabi Atlante reserves the right to make changes to the design and technical specifications of its products, with the authorization of ANVISA, in order to continue improving them.
 CONFIGURATION: CHAIR: NEW CROMA T5 • ANVISA: 10069210062 | WATER UNIT: SYNCRUS G2 • ANVISA 10069210063
 DELIVERY UNIT: Syncrus G3 / Syncrus G3 H • ANVISA 10069210075 | DENTAL LIGHT: Sirius G8 Sensor 3 LEDs • ANVISA 10069210069
 5100000174 - Version May 2023