

DENTAL UNIT

Prestige

DABI ATLANTE TECHNOLOGY **air**

DABI ATLANTE

- PNEUMATIC FOLDING ARMREST
- MULTIFUNCTION JOYSTICK PEDAL
- FLUSH SYSTEM
- WATER UNIT WITH PROXIMITY SENSOR AND CUP HOLDER
- 5 LED DENTAL LIGHT WITH PROXIMITY SENSOR (OPTIONAL)
- ANTI-CRUSHING SYSTEM
- RETRACTABLE RODS WITH LOCKING SYSTEM (HASTEFLEX VERSION)



ELECTRIC MICROMOTOR
*OPTIONAL



DENTAL UNIT

Prestige

hasteflex

DENTAL CHAIR

- Anti-crushing system
- Patient support arm with pneumatic activation
- Anatomic and bi-articulated headrest
- Joystick pedal with return to zero position, 3 working positions, spit position, and dental light activation.
- Trendelenburg position

DELIVERY UNIT

- Arm with pneumatic lock
- Retractable rods with locking system (Hasteflex Version)
- Removable stainless steel cover (Air Version)
- Flush system
- PAD panel with X-ray Viewer

WATER UNIT

- Easy access to the debris filters
- Water unit can be turned 60°
- Electric water activation
- Venturi and Vac Plus Suction Terminal
- Proximity Sensor
- Cup Holder

5 LED DENTAL LIGHT WITH PROXIMITY SENSOR

- Lower power consumption and longer durability
- Proximity sensor.
- Adjustable light intensity: 25.000 - 45.000 LUX
- Handles: Removable and autoclavable
- Frontal Protection: Removable

DABI COLOURS

UPHOLSTERY THAT INSPIRES YOUR WORK ROUTINE.

PVC OR VISCOELASTIC*



DARK GREEN



GOLD



GRAY



ORANGE



SILVER



NAVY BLUE



LIGHT GREEN



SKY BLUE



DESERT ROSE



NUT

LEATHER*



TERRACOTTA



NERO



MARINE



NUT



CARRARA

*Optional.

OPTIONS

Prestige

YOUR DENTAL OFFICE DESERVES EVOLUTION

N: NOT APPLICABLE | O: OPTIONAL | S: STANDARD

| DENTAL CHAIR | PRESTIGE AIR | PRESTIGE HASTEFLEX |
|---|--------------|--------------------|
| Column Lift System | N | N |
| Trendelenburg Position | S | S |
| Massager | O | O |
| Standard PVC | S | S |
| Viscoelastic PVC | O | O |
| Leather | O | O |
| Bi-articulated Headrest | S | S |
| Pneumatic Multi-articulated Headrest | N | N |
| 1x Fixed Armrest + 1x 90° Folding Armrest | N | N |
| 1x Fixed Armrest + 1x Pneumatic retractable Armrest | S | S |
| 1x Fixed Armrest + 1x Retractable Armrest (V Version) | N | N |
| Foldable Leg Support (at knee height) | N | N |
| Feet Protector | S | S |
| Foot Control (3 keys) | N | N |
| Foot Control (5 keys) | N | N |
| Foot Control (7 keys) | S | S |
| Foot Control Pedal - with Chip Blower function | N | N |
| Soft Comfort | O | O |
| Single Connection Box | O | O |
| Ambidextrous Dental Chair | N | N |

| DELIVERY UNIT | PRESTIGE AIR | PRESTIGE HASTEFLEX |
|--------------------------------------|-------------------|--------------------|
| Arm with Air Lock | S | S |
| Stainless Steel Top | S | N |
| Upper Tray (Acrylic) | N | N |
| Lower Tray (Injected) | N | S |
| Control Pad Kit | S | S |
| X-ray Viewer | S | S |
| Flush Kit | S | S |
| Manometer | S | S |
| High Speed Terminal | S | S |
| Low Rotation Terminal | N | N |
| Low Rotation Spray Terminal | S | S |
| Triple Syringe | S | S |
| Triple Injected Syringe | N | N |
| Syringe Heater Kit Equipment | O | O |
| Number of Terminals in the Equipment | UP TO 5 (3S & 2O) | UP TO 5 (3S & 2O) |
| Electric Micromotor Kit | O | O |
| Borden Terminal Kit - with Spray | O | O |
| Fibre Optic Terminal Kit | O | O |
| Curing Light Kit | O | O |
| Ultrasonic Scaler Kit | O | O |
| Bicarbonate Jet Kit | O | O |

| WATER UNIT | PRESTIGE AIR | PRESTIGE HASTEFLEX |
|---------------------------------------|--------------|--------------------|
| Proximity Sensor Kit | S | S |
| Cup Holder Kit | S | S |
| Electrical activation of water | S | S |
| Assistant Arm - without control panel | O | O |
| Assistant Arm - with control panel | O | O |
| Venturi Suction Terminal | S | S |
| Vac Plus Suction Terminal | S | S |
| Vacuum Pump Suction Terminal | O | O |
| Triple Injected Syringe | O | O |
| Triple Syringe Heating | N | N |
| Water Unit Syringe Heater Kit | O | O |

| DENTAL LIGHT | PRESTIGE AIR | PRESTIGE HASTEFLEX |
|----------------------------|--------------|--------------------|
| Persus LED | O | O |
| Sirius Sensor 3 LEDs | O | O |
| Sirius Sensor 5 LEDs | S | S |
| Sirius Sensor 5 (3x2) LEDs | O | O |

| OTHERS | PRESTIGE AIR | PRESTIGE HASTEFLEX |
|------------------------|--------------|--------------------|
| Multimedia Kit Support | O | O |

Illustrative Photos. Colours may have altered in the graphic reproduction. Dabi Atlante reserves the right to make changes to the design and technical specifications of its products, with the authorization of ANVISA, in order to continue improving them.
 CONFIGURATION: CHAIR: NEW CROMA • ANVISA: 10069210062 | WATER UNIT: SYNCRUS G3 • ANVISA 10069210063
 DELIVERY UNIT: Syncrus G3 / Syncrus G3 H • ANVISA 10069210075 | DENTAL LIGHT: Sirius G8 Sensor 5 LEDs • ANVISA 10069210069
 5100000174 - Version May 2023



dabiatlante



dabiatlante_oficial



dabiatlante

DABI ATLANTE

Dabi Atlante is a member of the Alliage group. | Contact: sales@alliage-international.com **alliage**