

OPTILIGHT MAX

CURING LIGHT

DABI/ATLANTE

MAX, is a wireless, high power and intensity LED curing light. With 3 pre-programmed light intensity modes, a beeper (in 5-second intervals), controls directly on the unit, and a practical and intelligent design - comfortable to operate. It's rechargeable battery and charging base ensures convenience and mobility.

FEATURES

- Wireless
- Transmission of light via 100% consistent fibre optic tip.
- High power and intensity LED.
- Operation controls on the device body.
- Allows adjustment of photoactivation times.
- 3 pre-programmed light intensity modes.
- Beep every 5 seconds.
- Rechargeable battery.
- Intelligent automatic shut off system.
- Resistant and easy to clean.
- Stand base with battery charger.

DENTIST BENEFITS

- Handles comfortably
- Movements with increased freedom, precision and lightness.
- No reduction in the quality of light intensity.
- Compact device - ergonomic and easy to install and transport.
- Can be used in various dental procedures, such as: restorative procedures, bracket bonding and activation of photoactivated materials such as sealants and base liners.

PATIENT BENEFITS

- No risk of dental heating, pulpal damage or discomfort for patients.
- Designed and built within the most advanced technology, to deliver results according to dental technical specifications.

INCLUDED ACCESSORIES:

Handpiece • Light guide • Eye protection • Charger base • Power supply



- 01 - Fibre Optic Tip
- 02 - Eyepiece Protector
- 03 - Control Panel
- 04 - Handpiece
- 05 - On/Off Button
(start/stop operation)
- 06 - Charging Base
- 07 - Display
- 08 - Time Adjustment Button
- 09 - Application Mode Selection
Button
- 10 - Power Supply
- 11 - Application Mode:
Continuous, Ramp & Pulsed

Illustrative photos.

Illustrative photos. Dabi Atlante reserves the right to make changes to the design and technical specifications of its products, with the authorization of ANVISA, in order to pursue its improvement. Anvisa: 10069210081. 51000000150 - Version May 2023



dabiatlante



dabiatlante_oficial



dabiatlante

DABI/ATLANTE

Dabi Atlante is a member of the Alliage group. | Contact: sales@alliage-international.com **alliage**