

SPECTRO 70X PERIAPICAL X-RAY

DABI ATLANTE



MOBILE COLUMN

- The Spectro 70X Periapical X-Ray - Mobile Column, has a base on 4 casters with high-stability columns.

WALL-MOUNTED WITH SCISSOR ARMS

- The Wall-Mounted Spectro 70X Periapical X-Ray, with scissor arms, allows greater reach and use in various positions.



- With arms that are lighter and easier to move, and a head with more precise positioning - even in the most challenging locations.
- Prepared to ensure patient and operator safety, with mechanisms that limit patient exposure.
- Tube (ampoule) with a 0.7 x 0.7mm focal point, aluminium filter >2.61mm and winding completely immersed in specialised oil.
- Rated voltage 70kVp and tube current 7.0 mA.

Illustrative Photos.



- 5 meter manual trigger - duly tested: Complying with the requirements of RDC No. 330, of December 20, 2019 and IN No. 95, of May 27, 2021 of ANVISA (Brazilian National Health Surveillance Agency).

Illustrative Photos. Dabi Atlante reserves the right to make changes to the design and technical specifications of the products, with the authorization of ANVISA, in order to pursue its improvement. Dental X-Ray Device: Anvisa 10069210087. 5100000249: Version May 2023