

CICLONE

1/2 HP and 1HP

DABI ATLANTE



SUCTION UNIT 1/2HP

Operates up to 2 dental sets

SUCTION UNIT 1HP

Operates up to 4 dental sets

Compact, with a high power suction system - for easy asepsis. Minimizes the possibility of cross-contamination and ensures cost savings.

FEATURES

- High power suction system.
- Aspires and drains all residuals from the oral area directly into the sewer, avoiding environmental contamination.
- Increases the biosecurity of the dental office.
- Thermal protection against motor temperature rise.
- Designed with corrosion-resistant materials, increasing the life-span of the components.
- Allows the regulation of suction power as needed (*BV suction kit - *Optional).
- Vacuum pump with baffle: reducing the noise level generated by the equipment.
- Made from bronze alloy.

BENEFITS

- Better performance: in prophylactic treatments and oral-maxillo surgeries.
- Time optimization
- Increased productivity:
 - for the user, due the patient needing to make use of the spittoon less frequently;
 - productivity for the clinic, by allowing the installation of two or four dental offices in simultaneous use.
- Infection control: Due to reduction of contaminated aerosols.
- Ease of operation: To operate the Vacuum Pump, simply remove the suction terminal from it's holder and it will automatically start operating. Return the suction terminal to its initial position, and the Vacuum Pump will turn off automatically.

Illustrative photos.

Illustrative photos. Dabi Atlante reserves the right to make changes in the design and technical specifications of products with the authorization of ANVISA in order to pursue its improvement. Anvisa: Non-correlated product. 51000000144 - Version May 2023



dabiatlante



dabiatlante_oficial



dabiatlante

DABI ATLANTE