

# VACUUM PUMP

## 1/2 HP and 1HP

**DABI**ATLANTE



### VACUUM PUMP ½HP

Allows to operate 2 dental sets

### VACUUM PUMP 1HP

Allows to operate 4 dental sets

The Vacuum Pump in the models ½ HP and 1HP has a high power suction system, besides being compact and of easy asepsis. The equipment minimizes the possibility of cross contamination in the office and ensures more savings for the dentist.

#### FEATURES

- High power suction system.
- Aspires and drains all residuals from the oral area directly into the sewer, avoiding environmental contamination.
- Increases the biosecurity of the dental office.
- It has thermal protection against motor temperature rise.
- Designed with corrosion-resistant materials, increasing the useful life of the components.
- Allows the regulation of the suction power as needed (\*BV suctor kit - \*Optional).
- Vacuum pump with baffle: reducing the noise level generated by the equipment.
- Made of bronze alloy.

#### DENTIST BENEFITS

- Better performance in prophylactic treatments and oral-maxillo surgeries.
- Optimization of work time.
- Increased productivity, because the patient goes less often to the spittoon.
- Infection control in the office, due to the reduction of contaminated aerosols.
- Increased productivity for the clinic, by allowing the installation of two or four dental offices in simultaneous use. rios en uso simultáneo.

To operate the Vacuum Pump, simply remove the suctor from the support and it will automatically start operating. When you put the suctor back to its initial position, the Vacuum Pump will turn off automatically.

Illustrative photos. Dabi Atlante reserves the right to make changes in the design and technical specifications of products with the authorization of ANVISA in order to pursue its improvement. Anvisa: Non-correlated product. 51000000144 - Version April/2021.

Illustrative photos.



dabiatlante



dabiatlante\_oficial



dabiatlante

**DABI**ATLANTE