

DENTAL SET

Versa

DABI ATLANTE TECHNOLOGY **air**

- MASSAGE SYSTEM
- COLUMN LIFT AND FOLDING LEG REST
- TRENDELENBURG POSITION
- MOBILE FOOT PEDAL WITH CHIP BLOWER FUNCTION
- FLUSH SYSTEM
- 5 LEDS DENTAL LIGHT WITH PROXIMITY SENSOR
- RETRACTABLE RODS WITH LOCKING SYSTEM (HASTEFLEX VERSION)

DABI ATLANTE



ELECTRIC MICROMOTOR (MME)
*OPTIONAL



DENTAL LIGHT 5 (3X2) LEDS*
*OPTIONAL

DENTAL SET

Versa hasteflex



DENTAL CHAIR

- Massage system.
- Column Lift and folding leg rest.
- Trendelenburg position.
- Folding armrest.
- Pneumatic multi-articulated headrest.
- Foot control with chip blower function, multifunctional activating chair, auxiliary unit, delivery unit and dental light.

DELIVERY UNIT

- Arm with pneumatic lock.
- Up to 6 terminals.
- Retractable rods with locking system (Hasteflex Version).
- Comes with 2 Trays (Air Version).
- Flush System.
- PAD panel with Xray Viewer

WATER UNIT

- Easy access to the debris filters.
- Electric activation of the water tank and cup holder with timer.
- Reach arm with PAD.
- Venturi and Vac Plus Suction Terminal.
- Proximity Sensor.
- Cupholder.

5 LEDS DENTAL LIGHT (STANDARD)

- Lower power consumption and longer durability.
- Proximity sensor.
- Selection of light intensity from 25,000 to 45,000 LUX.
- Handles: Removable and autoclavable.
- Frontal Protection: Removable

COLORS DABI

UPHOLSTERY THAT INSPIRES YOUR WORK ROUTINE.

VISCOELÁSTICO



DARK GREEN



GOLD



GRAY



ORANGE



SILVER



NAVY BLUE



LIGHT GREEN



SKY BLUE



DESERT ROSE



NUT

COURO*



TERRACOTTA



NERO



MARINE



NUT



CARRARA

*Optional.

OPTIONALS



YOUR DENTAL CHAIR DESERVES EVOLUTION

N: NOT APPLICABLE O: OPTIONAL | S: SERIES

DENTAL CHAIR	VERSA AIR	VERSA HASTEFLEX
Column Lift System	S	S
Trendelenburg Position	S	S
Massager	S	S
Standard PVC	N	N
Viscoelastic PVC	S	S
Leather	O	O
Bi-articulated headrest	N	N
Pneumatic Multi-articulated headrest	S	S
1 Fixed Arm Chair + 190° Folding Arm Chair	N	N
1 Fixed Arm Chair + 1 Pneumatic retractable Arm Chair	N	N
1 Fixed Arm Chair + 1 Retractable Arm Chair (V Version)	S	S
Foldable leg support at knee height	S	S
Feet Protector	S	S
Foot Control (3 keys)	N	N
Foot Control (5 keys)	N	N
Foot Control (7 keys)	N	N
Control panel with Chip Blower function (Mobile)	S	S
Soft Comfort	O	O
Single Connection Box	O	O
Ambidextrous Dental Chair	N	N

DENTAL UNIT	VERSA AIR	VERSA HASTEFLEX
Arm with Air Lock	S	S
Stainless steel top	N	N
Equipment Upper Tray (Acrylic)	S	N
Equipment Lower Tray (Injected)	N	S
PAD Kit (Electronic Controls)	S	S
X-ray viewer	S	S
Flush Kit	S	S
Manometer	S	S
High Speed Terminal	S	S
Low Rotation Terminal	N	N
Low Rotation Spray Terminal	S	S
Triple Syringe	N	N
Triple Injected Syringe	S	S
Syringe Heater Kit Equipment	O	O
Number of Terminals in the Equipment	Até 6 (3s e 3o)	Até 6 (3s e 3o)
Electric Micromotor Kit	O	O
Borden Terminal Kit with Spray	O	O
Optical Fiber Terminal Kit	O	O
Light Curing Kit	O	O
Ultrasound Kit	O	O
Bicarbonate Jet Kit	O	O

WATER UNIT	VERSA AIR	VERSA HASTEFLEX
Proximity Sensor Kit	S	S
Cup Holder Kit	S	S
Electrical activation of water	S	S
Swing arm without control panel	N	N
Swing arm with control panel	S	S
Venturi Sucker Terminal	S	S
Vac Plus Sucker Terminal	S	S
Vacuum pump suction terminal	O	O
Triple injected syringe	N	N
Triple Syringe Heating	O	O
UA Syringe Heater Kit	O	O

DENTAL LIGHT	VERSA AIR	VERSA HASTEFLEX
Persus L (Halogen)	N	N
Persus Led	N	N
Sirius Sensor 3 Leds	N	N
Sirius Sensor 5 Leds	S	S
Sirius Sensor 5 (3x2) Leds	O	O

OUTROS	VERSA AIR	VERSA HASTEFLEX
Multimedia Kit (Support and Monitor)	O	O

Illustrative Photos. Colors may have undergone alterations in the graphic reproduction. Dabi Atlante reserves the right to make changes in the design and technical specifications of its products, with the authorization of ANVISA, in order to continue improving them. CONFIGURATION: DENTAL CHAIR NEW VERSA: Anvisa 10069210062. SYNCRUS G3 WATER UNIT: ANVISA 10069210063. DELIVERY UNIT Syn Crus G8 FE and Syn Crus G8 HE: ANVISA 10069210075. DENTAL LIGHT Sirius G8 Sensor 5 Leds; Sirius G8 Sensor 5 (3x2) Leds: ANVISA 10069210069. S1000000174 - Version Jun/2021.



dabiatlante



dabiatlante_oficial



dabiatlante

DABI ATLANTE