

Discover the New Line
DENTAL UNIT 2024

CROMA LINE

- BIARTICULATED HEADREST
- 3 LED DENTAL LIGHT WITH PROXIMITY SENSOR OR PERSUS LED (CROMA 21)
- MULTIFUNCTION JOYSTICK
- INTEGRATED CONNECTION BOX
- ARTICULABLE DELIVERY UNIT ARM WITH PNEUMATIC LOCKING
- UP TO 5 TERMINALS ON THE EQUIPMENT, DEPENDING ON THE VERSION



CHAIR

- Single central articulation.
- Patient support arm foldable 90°.
- Bi-articulated headrest.
- Pedal with activation of the dental light, movements of the backrest and seat, back to zero and working positions programmable by the professional.
- Integrated connection box.
- Ambidextrous.

DELIVERY UNIT

- Arm with pneumatic lock
- Up to 5 terminals on the equipment, depending on the version
- Removable stainless steel top
- PAD panel that activates chair and dental light:
 - Seat movements
 - Backrest movements
 - Position back to zero
 - Working positions
 - Dental light

WATER UNIT

- Porcelain bowl, removable, with a drain to retain solids.
- Water unit foldable.
- Debris filter located in the base of the suction under the water unit body.
- Cup holder
- Proximity sensor
- Comes with: External and Cattani suction terminal

DENTAL LIGHT

- Sensor 3 LEDs (version Croma 21)
 - Proximity sensor.
 - Selection of light intensity from 10.000 to 30.000 LUX.
 - Handles: Removable and autoclavable.
 - Frontal Protection: Removable
- Persus LED (version Croma 11)
 - Monofocal
 - Multifaceted mirror
 - Intensity: 8.000 to 35.000 LUX (intensity choice by foot control)

SERIAL ITEMS



Persus Dental Light (LED light)

Monofocal for dental use with 1 LED optical system
 Multifaceted mirror with multicoating treatment
 Double protection of the mirror, in resistant, transparent material
 Bilateral handles in the shape of a handle which enable isolation, avoiding the risk of cross contamination
 Head made of resistant material, 620° swivel

Dental Light - 5 LEDs Sensor

- LED Lighting technology: Energy-efficient and long-lasting.
- Proximity sensor for activation and light intensity control.
- Cold light source, no heat generation.
- Rectangular spotlight for a larger work area.
- Innovative design with ergonomic handling.
- Rotatable head with durable material.
- Bilateral loop handles for isolation and autoclavable.
- Removable front protector for optical system protection.

- Sensor 3 LEDs
- Optical system with 3 LEDs
- Three Intensities: 10,000, 20,000 and 30,000 LUX

(only in version Cromia 2IM)



2 Fibre optic terminals

(only in version Cromia 2IM)



Control pad



Cup holder Proximity sensor

COMPARE THE VERSIONS




	CROMA 1IM	CROMA 2IM
Dental Chair	Bi-articulated headrest Folding armrest	Bi-articulated headrest Folding armrest
Delivery Unit	Air	Air
Delivery Unit's Accessories	Stainless steel tray Control pad	Stainless steel tray Control pad
Delivery Unit's Terminals	Triple syringe 2 High-speed terminals (Borden or Midwest) 1 Low-speed terminal (Borden or Midwest)	Triple syringe 1 High-speed terminals (Borden or Midwest) 1 Low-speed terminal (Borden or Midwest) 2 Fibre optic terminals
Water Unit's Accessories	Removable ceramic spittoon Cup holder Proximity sensor Selection valve	Removable ceramic spittoon Cup holder Proximity sensor Selection valve
Water Unit's Terminals	External suction terminal Cattani suction terminal	External suction terminal Cattani suction terminal
Dental Light	Persus LED	Sirius G8 3 LEDs

COLORS

DABI

Upholstery that inspires your work routine

PVC or Viscoelastic



Leather

