

Discover the New Line  
**DENTAL UNIT 2024**

**PERSONAL  
LINE**

BIARTICULATED HEADREST

FOLDING ARMREST

MULTIFUNCTIONAL JOYSTICK PEDAL

EQUIPMENT ARM WITH PNEUMATIC LOCK

WATER UNIT WITH PROXIMITY SENSOR  
AND CUP HOLDER

3 LED DENTAL LIGHT WITH PROXIMIT  
SENSOR

RETRACTABLE RODS WITH LOCKING SYSTEM  
(HASTEFLEX VERSION)



**CHAIR**

- Patient support arm foldable 90°.
- Anatomic and biarticulated headrest.
- Joystick pedal with back to zero position; 3 work positions; spit position, dental light activation and movement block.
- Ambidextrous.

**DELIVERY UNIT**

- Arm with pneumatic lock
- Up to 5 terminals on the equipment, depending on the version
- Removable stainless steel top
- PAD panel with negatoscope that activates chair and dental light:
  - Seat movements
  - Backrest movements
  - Position back to zero
  - Working positions
  - Dental light

**WATER UNIT**

- Debris filter located at the base of the chair.
- Water unit can be turned 90°.
- Electric activation of the water.
- External and Cattani suction terminal
- Proximity sensor.
- Cup holder.
- Assistant arm with Control pad

**DENTAL LIGHT**

- Sensor 3 LEDs
- Lower power consumption and longer durability.
- Proximity sensor.
- Selection of light intensity from 10.000 to 30.000 LUX.
- Handles: Removable and autoclavable.
- Frontal Protection: Removable

# SERIAL ITEMS

## Dental Light - 5 LEDs Sensor

- LED Lighting technology: Energy-efficient and long-lasting.
- Proximity sensor for activation and light intensity control.
- Cold light source, no heat generation.
- Rectangular spotlight for a larger work area.
- Innovative design with ergonomic handling.
- Rotatable head with durable material.
- Bilateral loop handles for isolation and autoclavable.
- Removable front protector for optical system protection.

- Sensor 3 LEDs
- Optical system with 3 LEDs
- Three Intensities: 10,000, 20,000 and 30,000 LUX



**2 Fibre optic terminals**  
except in versions Personal I and Personal IH



**Control pad with negatoscope**



**Assistant arm with Control pad**

## COMPARE THE VERSIONS




	PERSONAL I	PERSONAL IH	PERSONAL I2	PERSONAL I2H
<b>Dental Chair</b>	Bi-articulated headrest Folding armrest	Bi-articulated headrest Folding armrest	Bi-articulated headrest Folding armrest	Bi-articulated headrest Folding armrest
<b>Delivery Unit</b>	Air	Hasteflex	Air	Hasteflex
<b>Delivery Unit's Accessories</b>	Stainless steel tray Control pad with X-ray viewer Flush system	Lower tray Control pad with X-ray viewer Flush system	Stainless steel tray Control pad with X-ray viewer Flush system	Lower tray Control pad with X-ray viewer Flush system
<b>Delivery Unit's Terminals</b>	Triple syringe 2 High-speed terminals (Borden or Midwest) 1 Low-speed terminal (Borden or Midwest)	Triple syringe 2 High-speed terminals (Borden or Midwest) 1 Low-speed terminal (Borden or Midwest)	Triple syringe 1 High-speed terminals (Borden or Midwest) 1 Low-speed terminal (Borden or Midwest) 2 Fibre optic terminals	Triple syringe 1 High-speed terminals (Borden or Midwest) 1 Low-speed terminal (Borden or Midwest) 2 Fibre optic terminals
<b>Water Unit's Accessories</b>	Removable ceramic spittoon Assistant arm with Control pad Cup holder Proximity sensor Selection valve Master valve	Removable ceramic spittoon Assistant arm with Control pad Cup holder Proximity sensor Selection valve Master valve	Removable ceramic spittoon Assistant arm with Control pad Cup holder Proximity sensor Selection valve Master valve	Removable ceramic spittoon Assistant arm with Control pad Cup holder Proximity sensor Selection valve Master valve
<b>Water Unit's Terminals</b>	External suction terminal Cattani suction terminal	External suction terminal Cattani suction terminal	External suction terminal Cattani suction terminal	External suction terminal Cattani suction terminal
<b>Dental Light</b>	Sirius G8 3 LEDs	Sirius G8 3 LEDs	Sirius G8 3 LEDs	Sirius G8 3 LEDs

## COLORS

### DABI

Upholstery that inspires your work routine

#### PVC or Viscoelastic



#### Leather

