

Discover the New Line
DENTAL UNIT 2024

PRESTIGE
LINE

MULTIFUNCTION JOYSTICK PEDAL

FLUSH SYSTEM

WATER UNIT WITH PROXIMITY
SENSOR AND CUP HOLDER

5 LEDS DENTAL LIGHT WITH PROXIMITY
SENSOR

ANTI-CRUSHING SYSTEM

RETRACTABLE RODS WITH LOCKING SYSTEM
(PRESTIGE 11H AND PRESTIGE 12H)



CHAIR

- Anti-crushing system.
- Anatomic and biarticulated headrest.
- Joystick pedal with zero back position; 3 working positions; spit position and dental light activation
- Trendelenburg position.

DELIVERY UNIT

- Retractable rods with locking system (Hasteflex Version).
- Arm with pneumatic lock
- Up to 5 terminals on the equipment, depending on the version
- Removable stainless steel top
- PAD panel with negatoscope that activates chair and dental light

WATER UNIT

- Easy access to the debris filters.
- Water unit can be turned 60°.
- Assistant arm with Control pad
- Electric activation of the water tank and cup holder.
- External and Cattani Suction Terminal.
- Proximity sensor.
- Cup holder.

DENTAL LIGHT

- Sirius Sensor 5 LEDs
- Lower power consumption and longer durability.
- Proximity sensor.
- Selection of light intensity from 25.000 to 45.000 LUX.
- Handles: Removable and autoclavable.
- Frontal Protection: Removable

SERIAL ITEMS

Dental Light - 5 LEDs Sensor

- LED Lighting technology: Energy-efficient and long-lasting.
- Proximity sensor for activation and light intensity control.
- Cold light source, no heat generation.
- Rectangular spotlight for a larger work area.
- Innovative design with ergonomic handling.
- Bilateral loop handles for isolation and autoclavable.
- Removable front protector for optical system protection.

- Sensor 5 LEDs
- Optical system with 5 LEDs
- Three Intensities: 25,000, 35,000 and 45,000 LUX



Electric Micromotor

- Speed Transmission: 1:1
- Maximum Speed: 40,000 rpm
- Maximum Torque: 3.5 N.cm
- Irrigation System: Internal
- Lighting System: LED+
- Connection Type: ISO 3964 Short Coupling (INTRA)
- Weight: 57 grams

(except in versions Prestige I and Prestige IH)



Control pad with negatoscope

- Seat movements
- Backrest movements
- Blocking movements
- Position back to zero
- Activation of water from spittoon and cup holder
- Spit position and return to last position
- Working positions
- Dental light
- Pneumatic chair arm activation

Assistant arm with Control pad



COMPARE THE VERSIONS




	PRESTIGE I	PRESTIGE IH	PRESTIGE I2	PRESTIGE I2H
Dental Chair	Bi-articulated headrest Folding armrest	Bi-articulated headrest Folding armrest	Bi-articulated headrest Folding armrest	Bi-articulated headrest Folding armrest
Delivery Unit	Air	Hasteflex	Air	Hasteflex
Delivery Unit's Accessories	Stainless steel tray Control pad with X-ray viewer Flush system Manometer	Lower tray Control pad with X-ray viewer Flush system Manometer	Stainless steel tray Control pad with X-ray viewer Flush system Manometer	Lower tray Control pad with X-ray viewer Flush system Manometer
Delivery Unit's Terminals	Triple syringe 2 High-speed terminals (Borden or Midwest) 1 Low-speed terminal (Borden or Midwest)	Triple syringe 2 High-speed terminals (Borden or Midwest) 1 Low-speed terminal (Borden or Midwest)	Triple syringe 1 Low-speed terminal (Borden or Midwest) 2 Fibre optic terminals 1 Terminal with Electric Micromotor	Triple syringe 1 Low-speed terminal (Borden or Midwest) 2 Fibre optic terminals 1 Terminal with Electric Micromotor
Water Unit's Accessories	Removable ceramic spittoon Assistant arm with Control pad Cup holder Proximity sensor Selection valve Master valve	Removable ceramic spittoon Assistant arm with Control pad Cup holder Proximity sensor Selection valve Master valve	Removable ceramic spittoon Assistant arm with Control pad Cup holder Proximity sensor Selection valve Master valve	Removable ceramic spittoon Assistant arm with Control pad Cup holder Proximity sensor Selection valve Master valve
Water Unit's Terminals	External suction terminal Cattani suction terminal Triple syringe			
Dental Light	Sirius G8 5 LEDs			

COLORS

DABI

Upholstery that inspires your work routine

PVC or Viscoelastic



Leather

