

Discover the New Line
DENTAL UNIT 2024

VERSA
LINE

MASSAGE SYSTEM

COLUMN LIFT AND FOLDING LEGREST

TRENDELENBURG POSITION

MOBILE PEDAL WITH

CHIP BLOWER FUNCTION

FLUSH SYSTEM

3X2 LED REFLECTOR WITH

PROXIMITY SENSOR



CHAIR

- Massage system.
- Column Lift and leg support foldable.
- Trendelenburg position.
- Foldable support arm.
- Bi-articulated headrest
- Command pedal with chip blower function, multifunctional activating chair, auxiliary unit, equipment and reflector.

DELIVERY UNIT

- Arm with pneumatic lock.
- Up to 5 terminals on the equipment, depending on the version
- Retractable rods with system locking mechanism
(Versa 11H and Versa 12H).
- Removable stainless steel top
(versa 11 e versa 12).
- Flush System.
- PAD panel with negatoscope that activates chair and dental light
 - Seat movements
 - Backrest movements
 - Blocking movements
 - Position back to zero
 - Activation of water from spittoon and cup holder
 - Spit position and return to last position
 - Working positions

WATER UNIT

- Easy access to debris filters.
- Electric activation of the water in the tub and the coaster with timer.
- Arm reach with PAD.
- Terminal Sugador External and Cattani.
- Triple Syringe
- Proximity sensor.
- Cup holder.

DENTAL LIGHT

- Sirius 3x2 LED reflector
- Lower energy consumption and greater durability.
- Proximity sensor.
- Selection of luminous intensity from 15,000, 25,000 and 35,000 LUX (white light) 6,500 LUX (yellow light).
- Handles: Removable and autoclavable.
- Front Protector: Removable.

SERIAL ITEMS

Dental Light - 3x2 LED Sensor

- LED Lighting technology: Energy-efficient and long-lasting.
- Proximity sensor for activation and light intensity control.
- Cold light source, no heat generation.
- Rectangular spotlight for a larger work area.
- Innovative design with ergonomic handling.
- Bilateral loop handles for isolation and autoclavable.
- Removable front protector for optical system protection.

- Optical system with 5 LEDs: 3 white LEDs + 2 yellow LEDs
- Four Light Intensities – White Light: 15,000, 25,000 and 35,000 LUX and Yellow light: 6,500 LUX
- It allows the use of yellow light during clinical procedures with photoactivated materials, increasing their handling time by inhibiting the early polymerization of composite resins and other photoactivated materials



Electric Micromotor

- Speed Transmission: 1:1
- Maximum Speed: 40,000 rpm
- Maximum Torque: 3.5 N.cm
- Irrigation System: Internal
- Lighting System: LED+
- Connection Type: ISO 3964 Short Coupling (INTRA)
- Weight: 57 grams

(only in versaI2 and versa I2H)



Folding Leg Support

Massage System

COMPARE THE VERSIONS



	VERSA I	VERSA IH	VERSA I2	VERSA I2H
Dental Chair	Pneumatic headrest Folding armrest Massage System Column Lift Folding leg support	Pneumatic headrest Folding armrest Massage System Column Lift Folding leg support	Pneumatic headrest Folding armrest Massage System Column Lift Folding leg support	Pneumatic headrest Folding armrest Massage System Column Lift Folding leg support
Delivery Unit	Air	Hasteflex	Air	Hasteflex
Delivery Unit's Accessories	Stainless steel tray Control pad with X-ray viewer Flush system Manometer	Lower tray Control pad with X-ray viewer Flush system Manometer	Stainless steel tray Control pad with X-ray viewer Flush system Manometer	Lower tray Control pad with X-ray viewer Flush system Manometer
Delivery Unit's Terminals	Triple syringe 2 High-speed terminals (Borden or Midwest) 1 Low-speed terminal (Borden or Midwest)	Triple syringe 2 High-speed terminals (Borden or Midwest) 1 Low-speed terminal (Borden or Midwest)	Triple syringe 1 Low-speed terminal (Borden or Midwest) 2 Fibre optic terminals 1 Terminal with Electric Micromotor	Triple syringe 1 Low-speed terminal (Borden or Midwest) 2 Fibre optic terminals 1 Terminal with Electric Micromotor
Water Unit's Accessories	Removable ceramic spittoon Assistant arm with Control pad Cup holder Proximity sensor Selection valve Master valve	Removable ceramic spittoon Assistant arm with Control pad Cup holder Proximity sensor Selection valve Master valve	Removable ceramic spittoon Assistant arm with Control pad Cup holder Proximity sensor Selection valve Master valve	Removable ceramic spittoon Assistant arm with Control pad Cup holder Proximity sensor Selection valve Master valve
Water Unit's Terminals	External suction terminal Cattani suction terminal Triple syringe	External suction terminal Cattani suction terminal Triple syringe	External suction terminal Cattani suction terminal Triple syringe	External suction terminal Cattani suction terminal Triple syringe
Dental Light	Sirius G8 3x2 LEDs	Sirius G8 3x2 LEDs	Sirius G8 3x2 LEDs	Sirius G8 3x2 LEDs

COLORS

DABI

Upholstery that inspires your work routine

PVC or Viscoelastic



Leather

